

# Spray application of *in situ* cross-linkable multi-layered hydrogels for medical applications

---

## 医療用*in situ*架橋マルチレイヤーハイドロゲルの スプレー投与システム


医科学専攻修士課程修士論文発表会

指導教員：伊藤大知 准教授

41-156214 Zheng Ying Grace

# In situ cross-linkable hydrogels

- Gels formed by simple mixing of two different polymers
- Merits:
  - ✓ Cover complex tissue geometries
  - ✓ Used in both open & laparoscopic surgeries
- Used in various biomedical applications:
  1. Scaffolds for tissue engineering
  2. Clinical procedures
    - Hemostasis
    - Tissue adhesion prevention



# Biomedical spray devices for clinical uses

- One method of delivery of *in situ* cross-linkable hydrogels
- Advantages:
  - ✓ Even delivery onto varying tissue surfaces
  - ✓ Less tissue handling

## Fibrin sealants

Bolheal®




Tisseel®



>50% of hemostats market share

## Tissue adhesion barriers

Adspray®




Few

**Sprays for adhesion prevention remains largely unstudied**

# Tissue adhesion

- Abnormal fibrous bands attaching 2 adjacent body tissues
- Occurs after 54% of abdominal surgeries
- Problems such as pelvic pain etc.



G.S.Dizerega and J.D. Campeau: Human Reproduction  
Update 2001(7)547-555

## Tissue adhesion barriers

- Physical separation of tissues
- Commercially available adhesion barriers have limitations:

**Seprafilm®**



**Interceed®**



**X Brittle**

**X Low efficacy in surgeries with heavy bleeding**

**Limited efficacy of a single material**


<http://www.colmed.in/interceed-adhesion-prevention/4350-interceed-adhesion-barrier-3-in-x-4-in-box-of-10.html>

<https://www.classactionlawsuitshelp.com/seprafilm-class-action-lawsuit/>

**Need for improved material → Combination of various materials**


# Objective

Spray application of multi-layered *in situ* cross-linkable hydrogels for tissue adhesion prevention



# Materials

Viscosity of materials:



	Liquid A	Liquid B
<b>Top</b> (away from wound surface)	1wt% Alginate	3wt% Carbomethyl cellulose (CMC)+ 0.1M CaCl <sub>2</sub> <b>Alg-CMC</b>
<b>Bottom</b> (in contact with wound surface)	2wt% Alginate+ 1wt% Transglutaminase (TG)*	20wt% Gelatin+0.1M CaCl <sub>2</sub> <b>Alg-gelatin</b>

\* [TG] was evaluated by fellow lab member Kikuchi E. 東京大学工学部卒論 (2016)

# Experimental flow

## Spray characterization

- Spray coverage
- Gel thickness
- Droplet size

**Optimization of spray conditions**

## *In vitro* analysis

- Cell viability
- Cell migration

**Effect of spray conditions and materials on cells**



## *In vivo* analysis

- Rat hepectatotomy adhesion model


**Response to multi-layered hydrogels**

# Spray characterization

**N<sub>2</sub> gas flow rate:  
3, 5, 7L/min**




Spray coverage  
Gel thickness



**Homogeneity**

Droplet sizes




While stirring

**Mixing efficiency**



# In vitro evaluation of cell damage by spraying


Cell types: NIH-3T3 fibroblasts, MeT-5A mesothelial cells  
Gas flow rate: 3, 5, 7L/min  
Spray distance: 5, 10, 20cm




# *In vitro* evaluation of cell migration

- To functionally differentiate the two layers

## 'Scratch' wound healing assay





## Cell penetration assay



# In vivo evaluation


## ➤ Rat hepatectomy model

- 8-11 week old SD rat ♂
- **3cm** of left lateral lobe of liver was transected




Evaluation after 1 week

### EXTENT




- Control
- Alg-gelatin
- Alg-CMC
- Double layer




### GRADE

- Grade 0: No adhesion
- Grade 1: Spontaneously separating adhesion
- Grade 2: Adhesions separated by blunt dissection
- Grade 3: Adhesions separated by sharp dissection

# Evaluation of homogeneity



## Evaluation of homogeneity




- Coefficient of variation (COV) =  $\frac{\text{Standard deviation } (\sigma)}{\text{Mean } (\mu)}$  of gel thickness


**The higher the gas flow rate and spray distance, the higher the homogeneity**

# Droplet size measurements


Spray distance: 5cm  
Gas flow rate: 3L/min



Spray distance: 10cm  
Gas flow rate: 7L/min



Material: 0.5, 1, 2, 3wt% Alginate



Smaller droplets with larger spray coverage

**High gas flow rates**  
**High spray distances**





**High homogeneity**  
**High mixing efficiency of gels**




# In vitro analysis: WST assay

## ➤ Alg-CMC at varying spray parameters




Materials sprayed at **high gas flow rate** and **short distances** caused significant cell death



# In vitro analysis: Live/Dead Staining

Spray distance = 5cm

Spray distance = 20cm




Materials sprayed at high gas flow rate and short distances caused **cell detachment** and **cell death of remaining attached cells**





# In vitro analysis: 'scratch' assay

## ➤ Observations of 'scratch' wound closure




**Alg-gelatin promoted horizontal cell migration**





# Fabrication of multi-layered hydrogels

- Different colored food dyes added to materials before spraying
- Transverse section of gel:




**Multi-layered hydrogel was successfully fabricated**

# *In vitro* analysis: cell penetration assay


➤ Transverse gel sections after 24h

➤ Fluorescence intensity


Migration  
direction




Alg-CMC layer




Migration  
direction




Alg-gelatin layer



Migration  
direction




Double layer



Alg-  
gelatin

Alg-CMC




**Cells migrated quickly  
through Alg-gelatin,  
not Alg-CMC**


# In vivo analysis

➤ Cut liver surface-greater omentum (離断面と大網)

Adhesion grade




Adhesion extent



- Alg-gelatin single layer worsened adhesion extent
- Double layer improved adhesion extent
- No significant difference between control and double layer


# Discussion

- At LOW gas flow rate:




- At HIGH gas flow rate:

Droplet size ↓  
Spray coverage ↑



- ✓ Reduced adhesion formation



- ✓ Promoted wound healing

**X Worsened adhesion formation**

Mixing efficiency  
Homogeneity

**Trade-off**

Cellular damage

**Precise control of thickness  
of each layer is needed**

# Conclusions

- Spray distribution and droplet sizes were measured
- Impact of spray conditions on cell viability was evaluated
- The wound healing property, and the anti-adhesive property of the respective hydrogel layers were shown
- Multi-layered hydrogel was successfully fabricated
- Efficacy of the multi-layered hydrogel delivered by spraying was evaluated using the rat hepatectomy model

# Acknowledgements

Special thanks to Kaketsuken/Sumitomo Bakelite for the provision of double syringes, Kimica for the provision of alginate, Nitta Gelatin for the provision of gelatin, and Ajinomoto for the provision of transglutaminase.

Prof. Kokudo, Associate Prof. Hasegawa, Dr. Inagaki, Dr. Toda from the HepatoBiliary-Pancreatic Surgery Division, Department of Surgery, Graduate school of Medicine at the University of Tokyo.

Associate Prof. Ito for all the guidance, Assistant Prof. Ohta for help with the animal experiments, and fellow lab members for their advice and support in one way or another. 😊